


ELISA Buffer Kit & ELISA Development Kits (EDK)

Each ELISA Development Kit (EDK) contains the following components (RED):
 Capture antibody, Detection antibody (biotinylated), a calibrated Standard and Avidin-HRP Conjugate

Standard EDK = sufficient for 1000 ELISA plate wells
Mini EDK = sufficient 200 ELISA plate wells

The other necessary components (GRAY) to assay ELISA can be purchased as **ABTS ELISA Buffer Kit** (cat#: 900-K00) or **TMB ELISA Buffer Kit** (cat#: 900-T00)



EDK Kit Contents:


- Standard
- Capture Antibody
- Detection Antibody
- Avidin-HRP Conjugate

EDK + ELISA Buffer kit = Complete ELISA Kit

ELISA Buffer Kit:

- 10 Uncoated Microplates
- 20x Plate Coating Buffer
- 20x Wash Buffer
- 1x Blocking Buffer
- 20x Sample Diluent
- Ready-to-use Substrate (ABTS or TMB)
- 50 Plate Sealing Films

*The ELISA Buffer kit contains all of the necessary components (Gray) to assay 10 ELISA plates.



Search Your Protein here

<https://www.peprotech.com/en-US/Pages/ProductIndex/Detection-Products/Standard-ELISA-Development-Kits>

